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**OLM INTRACRANIAL  
PRESSURE  
MONITORING KIT  
MODEL 110-4B**

---

**DIRECTIONS FOR USE**

Read carefully prior to use

---

**WARNING:**

This product contains a chemical (Ethylene Oxide) known to the State of California to cause cancer, birth defects or other reproductive harm.

**CAUTION: FEDERAL (U.S.A.) LAW RESTRICTS THIS DEVICE  
TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

 **INTEGRA™**  
Neuro Sciences  
**CAMINO®**


# OLM\* INTRACRANIAL PRESSURE MONITORING KIT

## MODEL 110-4B

\* The OLM ICP Kit was developed in cooperation with Richard C. Ostrup, M.D., Thomas G. Luerssen, M.D. and Lawrence F. Marshall, M.D. (San Diego, California)

ENGLISH

Figure 1

## SYSTEM DESCRIPTION

The OLM Intracranial Pressure Monitoring Kit consists of a sterile transducer-tipped pressure monitoring catheter with accessory items to be used as a diagnostic tool for rapidly determining and continuously monitoring intracranial pressure. Since its method of measuring pressure is unique, please read this section carefully.

Unlike ordinary pressure monitoring systems, the Camino® Catheter has a miniature transducer at the distal tip. This unique design eliminates the need for a "fluid-filled system" to carry pressure waves to an external transducer. Thus, the problems inherent in such fluid-filled systems are eliminated. The result is a precise pressure measurement and an artifact free, high fidelity waveform trace.

The OLM Intracranial Pressure Monitoring Kit contains the following accessory items for use with the Camino Catheter (Figure 1);

- Camino Bolt with Compression Cap and Turning Wings,
- .106 inches (2.71 mm) diameter No. 36 Drill Bit with Safety Stop,
- Spacer to adjust seating depth of the bolt,
- Strain Relief—Protective Sheath,
- .050 inches (1.3 mm) diameter Stylet,
- Hex wrench,
- Zero adjustment tool

Note: A hand drill and various standard surgical instruments and supplies are required to place the ICP monitoring system.

A set of instruments and supplies is available from Integra NeuroSciences in our Cranial

Access Kits.

## INDICATIONS

The use of the OLM\* Intracranial Pressure Monitoring Kit by a qualified neurosurgeon is indicated when direct measurement of intracranial pressure in the parenchyma or the subarachnoid space, is clinically important.

## CONTRAINDICATIONS

This device is not intended for any use other than that indicated.

## RISK AND COMPLICATIONS

Maintenance of sterility during placement and subsequent handling of the Camino Bolt is essential. Sterile technique should be used at all times when inserting, correcting, or adjusting the Camino Bolt and the Camino Transducer-Tipped Pressure Monitoring Catheter.

If monitoring is continued for more than 5 days, placement of a new system under sterile conditions is recommended.

Placement of the Camino Bolt must be carried out by a qualified neurosurgeon to avoid penetration of the cerebral cortical surface during the drilling process, which requires standard surgical procedure and skill.

Hemorrhage from the dura or cortical surface at the bolt insertion site may occur. Patients should be tested for normal blood clotting function prior to bolt placement. Determination of possible extracerebral, subarachnoid, or intracerebral hemorrhage at the bolt insertion site will be the

responsibility solely of the operating neurosurgeon. Appropriate steps and proceedings to control such hemorrhage should be taken when indicated by the neurosurgeon in charge.

## PRECAUTIONS


- Extreme bending and/or kinks can impair the performance of the Fiber Optic Pressure Transducer. Exercise caution when handling the catheter.
- The catheter is designed for **SINGLE USE ONLY. DO NOT RESTERILIZE OR REUSE.** Camino 110-4B is supplied in a single use package and is guaranteed to be sterile and non-pyrogenic unless opened or damaged. The product is intended for use as an intracranial pressure monitor and is not to be reused. Any attempts to resterilize or reuse the product will damage the catheter and impair its ability to function as intended.
- Use aseptic technique throughout procedures.
- Maintain the insertion site with regular meticulous redressing using aseptic technique.
- Do not attach anything to transducer air vent. Vent must remain open for proper operation (Figure 1).

## CAUTION

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

## INSTRUCTIONS FOR USE

### INSERTION METHOD FOR THE CAMINO ICP BOLT

- Area of insertion: The standard right and left prefrontal areas are the primary areas of insertion. This region allows the patient to have his head rotated from side to side and still remain in a supine position without interference with the intracranial pressure monitoring function. In addition, the incision will be carried behind the hairline in the majority of patients and therefore be cosmetically acceptable (Figures A and B).
- After the insertion site has been chosen, the area is shaved and prepped in a sterile

fashion, usually with a betadine solution. The shaved and prepared area is then draped with sterile towels. The area of the incision, which usually lies two to three centimeters anterior to the coronal suture in the mid-pupillary line is infiltrated subcutaneously with 1% Xylocaine. An approximately half centimeter linear incision is made and carried to the bone. A small mastoid type of retractor is then inserted to provide a good bone exposure and hemostasis of the skin edges.

- The safety stop on the drill bit provided in the kit can be positioned as desired by loosening the setscrew with the hex wrench, sliding the stop into position, and retightening the setscrew.
- The drill bit is then secured to a twist drill and, in a standard fashion, a twist drill hole is made through the outer and inner tables of the skull. The surgeon needs to be careful when penetrating the inner table to minimize any potential for parenchymal injury.
- After penetration of the inner table, the drill is removed and the hole is irrigated with sterile saline. An 18G spinal needle is then used to open the dura in a cruciate fashion. The stylet can be inserted to ensure adequate opening of the dura.
- Following opening of the dura, the Camino Bolt is screwed manually into the skull. The seating depth of the Camino Bolt will be at the surgeon's discretion pending the thickness of the skull. This will be approximately 2-3 mm for the neonatal age group, 3-5 mm for the pediatric age group and 5 mm to 1 cm for adults. If desired, the spacer can be used as a guide, otherwise the spacer can be removed and discarded.
- The stylet provided in the kit is inserted through the Camino Bolt and the dura to clear the passage for the Camino Transducer-Tipped Catheter.
- The Camino Bolt is irrigated with non-bacteriostatic sterile saline.

### CAMINO TRANSDUCER-TIPPED PRESSURE MONITORING CATHETER PREPARATION PRIOR TO INSERTION

- The Camino Catheter can be used with the Camino V420/420 series or with Camino MPM series pressure monitors. For Camino monitor set-up and use, refer to the appropriate Camino monitor Directions For Use.
- Remove the Camino Catheter from its sterile package and firmly attach the transducer connector to the preamp connector (Figure 2). If the monitor display does not read zero after a short system self-check delay, use the tool from the catheter kit to turn the zero adjustment on the bottom

side of the transducer connector until the monitor display reads zero (Figure 3).

Figure 2


Figure 3



Figure 4

- **WHEN THE RED DEPTH INDICATOR IS VISIBLE ABOVE THE STRAIN RELIEF, CATHETER TIP (TRANSDUCER) POSITION MUST BE CHECKED TO ENSURE DESIRED DEPTH IS ACHIEVED.**

- The surgeon may easily vary the insertion depth by locating his fingers at the proper cm mark before performing the above step. For example, placing the fingers at 5.5 cm will locate the tip of the catheter 1 cm beyond the end of the bolt, into the parenchyma.

**IT IS RECOMMENDED THAT THE CATHETER BE DISCONNECTED FROM THE PRE-AMP CONNECTOR WHEN THE PATIENT IS MOVED. THIS WILL NOT AFFECT CALIBRATION.**

## CONTINUOUS PRESSURE MONITORING


Since the Camino Catheter has a miniaturized transducer at the distal tip, it requires no fluid-filled system. Thus, the need for an external transducer, pressure dome and pressure tubing is eliminated. As a result, pressure may be monitored continuously without recalibration.

## INSERTION OF CAMINO CATHETER


- To measure intracranial pressure, insert the Camino Catheter into the Camino Bolt, using the cm markings on the catheter to gauge insertion depth (Figure 4). If the surgeon places his fingers at the 5 cm mark (double dot), then inserts the catheter until his fingers touch the top of the bolt, the tip of the catheter will be 0.5 cm beyond the end of the bolt, approximately in the subarachnoid space. He should pull the catheter back slightly, then turn the compression cap on the Camino Bolt clockwise to secure the catheter in place. If using a bedside monitor and/or Camino monitor, verify pressure waveform. If necessary, loosen the compression cap, reposition the Camino Catheter, and retighten the compression cap.
- While holding the catheter from above, slide the strain relief sheath down and secure it onto the compression cap.

# RECOMMENDED INSERTION PROCEDURE SUMMARY


After drilling hole in skull and penetrating the dura,


1. Screw in bolt.


2. Clear path with stylet.


3. Irrigate.


4. **Zero the transducer**, then place fingers at 5 cm and insert catheter until fingers touch top of bolt.

5. Pull back slightly.


6. Secure by tightening compression cap. Verify waveform, reposition if necessary.

7. Hold catheter from above, pull strain relief down and attach to bolt.

7a. When the red depth indicator is visible above the strain relief, catheter tip (transducer) position must be checked to ensure desired depth is achieved.

8. To **protect fiber optics**, position as shown and secure preamp connector with clip.

## INTERNATIONAL MARKINGS

### CATHETER EXTENSION CABLE CONNECTIONS

Symbol also signifies:

**TYPE BF EQUIPMENT:** Protected against electric shock, having an F-type isolated (floating) applied part.

**STERILE EO** Sterilized using ethylene oxide.

Date of expiration.

Do not reuse, single use device.

Manufacturer.

Do not resterilize.

**TYPE CF EQUIPMENT:** Provided with a degree of protection against electric shock higher than Type BF equipment, having an F-type isolated (floating) applied part; intended for direct cardiac application.

Refer to Directions for Use.

Lot Number.

**Rx ONLY**

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.**

**Authorized representative in the European Community.**

## SPECIFICATIONS

Transducer size	4F (1.35mm Diameter)
Transducer type	Fiber optic
Frequency response (system)	
Model MPM-1/SPM-1	100Hz (-3dB)
Model V420/M420/420XP	120Hz (-3dB)
Model 420	33Hz (-3dB)
Measurement range (system)	-10 to 125 mmHg
Zero drift (system)	
First 24 hours (maximum)	0 ± 2 mmHg
5 days (typical)	less than ± 1 mmHg per day
Temperature coefficient	Max of 3 mmHg over temperature range of 22°C to 38°C (70°F-100°F)
Linearity and hysteresis (system)	
Pressure range:	
-10 to 50 mmHg	± 2 mmHg or better
51 to 125 mmHg	± 6% of reading or better
Reference pressure	Atmosphere
Overpressure	-700 to 1250 mmHg

Camino and Integra NeuroSciences are registered trademarks at Integra LifeSciences Corporation in the United States and/or other countries. The Integra wave logo is a trademark of Integra LifeSciences Corporation © 2010 Integra LifeSciences Corporation. All Rights Reserved.

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Integra NeuroSciences Ltd.,